

PACKAGES AND ACCESSORIES

1. Plastic bottles for liquids	864
2. Plastic canisters for liquids	865
3. Plastic bottles for solids	865
4. Glas bottles for liquids	865
5. Glas bottles for solids	866
6. Transport packagings / fall protection packagings	866
6.1 Polystyrene boxes	866
6.2 Hazardous goods cartons	867
7. Accessories	867

1. Plastic bottles for liquids

0.1 Plastic bottles, natural-colored	product code	min. quantity packaging unit	1× EUR
with DIN 32 lock	50000.0000	10 piece	1,70
0.25 Plastic bottles, natural-colored	product code	min. quantity packaging unit	1× EUR
with DIN 32 lock	50000.1000	10 piece	2,00
0.5 Plastic bottles, natural-colored	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.2000	4 piece	3,60
0.5 Plastic bottles, natural-colored - degassing	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 degassing lock	50000.2010	4 piece	5,50
0.5 Plastic bottles, black	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.2100	4 piece	3,90
1 Plastic bottles, natural-colored	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.3000	4 piece	4,50
1 Plastic bottles, natural-colored - degassing	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 degassing lock	50000.3010	4 piece	6,50
1 Plastic bottles, black	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.3100	4 piece	5,10
2.5 Plastic bottles, natural-colored	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.4000	4 piece	8,20
2.5 Plastic bottles, natural-colored - degassing	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 degassing lock	50000.4010	4 piece	10,00
2.5 Plastic bottles, black	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.4100	4 piece	9,60
5 Plastic bottles, natural-colored	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.5300	4 piece	12,30
5 Plastic bottles, natural-colored - degassing	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 degassing lock	50000.5310	4 piece	13,40
5 Plastic bottles black	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 45 lock	50000.5100	4 piece	13,80

2. Plastic canisters for liquids

5 l Plastic canisters, natural-colored	product code	min. quantity packaging unit	1× EUR
UN-approved with lock no. 51	50000.5020	4 piece	8,20

5 l Plastic canisters, natural-colored - degassing	product code	min. quantity packaging unit	1× EUR
UN-approved with degassing lock no. 51	50000.5010	4 piece	9,70

5 l Plastic canisters, blue	product code	min. quantity packaging unit	1× EUR
UN-approved with lock no. 51	50000.5120	4 piece	8,20

10 l Plastic canisters, natural-colored	product code	min. quantity packaging unit	1× EUR
UN-approved with lock no. 51	50000.6020	1 piece	10,60

10 l Plastic canisters, natural-colored - degassing	product code	min. quantity packaging unit	1× EUR
UN-approved with degassing lock no. 51	50000.6010	1 piece	12,10

10 l Plastic canisters, blue	product code	min. quantity packaging unit	1× EUR
UN-approved with lock no. 51	50000.6120	1 piece	10,60

3. Plastic bottles for solids

0.25 l Plastic bottles wide neck, white	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 60 lock	50000.1600	10 piece	2,30

0.75 l Plastic bottles wide neck, white	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 60 lock	50000.2600	4 piece	3,30

1.5 l Plastic bottles wide neck, white	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 80 lock	50000.3600	4 piece	4,60

2.5 l Plastic bottles wide neck, white	product code	min. quantity packaging unit	1× EUR
UN-approved with DIN 80 lock	50000.4600	4 piece	5,70

4. Glas bottles for liquids

0.1 l Glas bottles	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 32 lock	50000.0700	10 piece	2,40

0.25 l Glas bottles	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 32 lock	50000.1700	10 piece	3,10

0.5 l Glas bottles	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 45 lock	50000.2700	4 piece	5,40
1 l Glas bottles	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 45 lock	50000.3700	4 piece	6,70
2.5 l Glas bottles	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 45 lock	50000.4700	4 piece	12,20
2.5 l Glas bottles with Surlyn®-splinter protection-coating	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 45 lock	50000.4750	4 piece	48,50

5. Glas bottles for solids

0.25 l Glas bottles wide neck	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 45 W lock	50000.1800	10 piece	2,30
0.5 l Glas bottles wide neck	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 54 W lock	50000.2800	4 piece	4,90
1 l Glas bottles wide neck	product code	min. quantity packaging unit	1× EUR
brown glass with DIN 54 W lock	50000.3800	4 piece	6,10

6. Transport packagings / fall protection packagings

6.1 Polystyrene boxes

Polystyrene box for 1 x 0.1 l glas bottle	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.0711	10 piece	1,20
Polystyrene box for 6 x 0.1 l glas bottles	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.0716	1 piece	4,10
Polystyrene box for 1 x 0.25 l glas bottle	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.1711	10 piece	1,60
Polystyrene box for 6 x 0.25 l glas bottles	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.1716	1 piece	5,80
Polystyrene box for 1 x 0.5 l glas bottle	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.2711	4 piece	2,50
Polystyrene box for 6 x 0.5 l glas bottles	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.2716	1 piece	8,90
Polystyrene box for 1 x 2.5 l glas bottle	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.4711	1 piece	7,40

Polystyrene box for 1 x 1 l glas bottle	product code	min. quantity packaging unit	1× EUR
UN-approved	51000.3711	4 piece	3,10

6.2 Hazardous goods cartons

Hazardous goods carton for 1 x 1 plastic bottle	product code	min. quantity packaging unit	1× EUR
corrugated folding carton, UN-approved suitable for 1 plastic container containing up to 1 l	51000.3001	10 piece	2,20

Hazardous goods carton suitable for 4 x 1 plastic bottles	product code	min. quantity packaging unit	1× EUR
corrugated folding carton, UN-approved suitable for 4 plastic containers containing up to 1 l	51000.3004	10 piece	3,10

Hazardous goods carton suitable for 12 x 1 plastic bottles	product code	min. quantity packaging unit	1× EUR
corrugated folding carton, UN-approved suitable for 12 plastic containers containing up to 1 l	51000.3012	10 piece	5,10

Hazardous goods carton for 4 x 1 glas bottles	product code	min. quantity packaging unit	1× EUR
corrugated folding carton, UN-approved suitable for 4 x 1 l glas bottles for liquids	51000.3704	1 piece	7,50

Hazardous goods carton suitable for 4 x 2.5 plastic bottles	product code	min. quantity packaging unit	1× EUR
corrugated folding carton, authorized by UN suitable for 4 plastic containers containing up to 2.5 l	51000.4004	10 piece	3,50

Hazardous goods carton for 4 x 2.5 glas bottles	product code	min. quantity packaging unit	1× EUR
corrugated folding carton, UN-approved suitable for 4 x 2.5 l glas bottles for liquids	51000.4704	1 piece	10,00

7. Accessories

Hand pumping ball	product code	min. quantity packaging unit	1× EUR
for the withdrawal system 06280 for 25 l Ekm for solvents with PTFE sealings	06281.0900	1 piece	74,00

Key for IBC cap DN 150	product code	min. quantity packaging unit	1× EUR
Key for IBC cap DN 150	25466.0900	1 piece	76,50

Key for Labcube® cap	product code	min. quantity packaging unit	1× EUR
Key for Labcube® cap	25465.0900	1 piece	3,90

Key for bottle locks	product code	min. quantity packaging unit	1× EUR
for caps according to DIN	06276.0900	1 piece	#NV

Key for Pka and Pkb 5 l / 10 l	product code	min. quantity packaging unit	1× EUR
for canister locks no. 51	06272.0900	1 piece	3,90

Key for Pka and Pkb 20 l / 30 l	product code	min. quantity packaging unit	1× EUR
for canister locks no. 61	06273.0900	1 piece	5,70

Key for 25 l Blk and other barrels	product code	min. quantity packaging unit	1× EUR
bronze version, sparkling excluded	06278.0900	1 piece	108,50
Key for 60 l Pka	product code	min. quantity packaging unit	1× EUR
for canister locks no. 71	12850.0900	1 piece	12,50
Kraft spout for Pka and Pkb 5 l / 10 l	product code	min. quantity packaging unit	1× EUR
spout for canister locks no. 51	06274.0900	1 piece	6,20
Kraft spout for Pka and Pkb 20 l / 30 l	product code	min. quantity packaging unit	1× EUR
spout for canister locks no. 61	06275.0900	1 piece	6,60
Kraft spout for Pka and Pkb 60 l	product code	min. quantity packaging unit	1× EUR
spout for canister locks no. 71	28451.0900	1 piece	7,20
Kraft spout brass	product code	min. quantity packaging unit	1× EUR
brass spout for screwing in for 30 l iron pots with 3/4" connection threads (the spout is self-closing)	06279.0900	1 piece	30,50
Kraft take-off 5 l / 10 l	product code	min. quantity packaging unit	1× EUR
for the transition of 5 l and 10 l canisters (Pka/Pkb) on GL 45	15979.0900	1 piece	17,30
Kraft take-off 20 l / 25 l / 30 l	product code	min. quantity packaging unit	1× EUR
for the transition of 20 l, 25 l and 30 l canisters (Pka/Pkb) on GL 45	17384.0900	1 piece	17,30
Kraft take-off Labcube®	product code	min. quantity packaging unit	1× EUR
connection adapter for Labcube® for connecting with titrators by Mettler	50050	1 piece	6,60
Spouts 5 l / 10 l	product code	min. quantity packaging unit	1× EUR
for Pka / Pkb canister locks no. 51	06270.0900	1 piece	3,10
Spouts 20 l / 30 l	product code	min. quantity packaging unit	1× EUR
for Pka / Pkb canister locks no. 61	06271.0900	1 piece	3,10
Standkraft	product code	min. quantity packaging unit	1× EUR
tapping tablet for Labcubes®	26798.0900	1 piece	81,50
Withdrawal system – stainless steel	product code	min. quantity packaging unit	1× EUR
for 25 l Ekm for solvents with PTFE sealings	06280.0900	1 piece	503,00