

CONTAINER ABBREVIATIONS

Flexibility is one of our greatest strengths. This is reflected consistently in our range of packaging solutions to meet customer requirements.

The following list gives all the package keys in our portfolio, their abbreviations and the corresponding type of packaging.

If you cannot find the packaging you need in our catalogue, please don't hesitate to contact us. We would be pleased to make you an individual offer.

Package key	Package size	Package unit	Package article	Package article in full	Thread
.0000	0.1	l	Pfl	Plastic bottle	DIN 32
.0300	0.1	l	PFA	Perfluorokoxy bottle - reusable	23 mm *
.0600	0.1	kg	Pfw	Plastic bottle, wide neck	DIN 60
.0700	0.1	l	Fl	Glass bottle	DIN 32
.0800	0.1	kg	Gl	Glass bottle, wide neck	DIN 40
.1000	0.25	l	Pfl	Plastic bottle	DIN 32
.1060	0.25	l	TNB	Twin-neck bottle	24 mm *
.1300	0.25	l	PFA	Perfluorokoxy bottle - reusable	23 mm *
.1600	0.25	kg	Pfw	Plastic bottle, wide neck	DIN 60
.1700	0.25	l	Fl	Glass bottle	DIN 32
.1800	0.25	kg	Gl	Glass bottle, wide neck	DIN 45W
.2000	0.5	l	Sfl	Safety plastic bottle	DIN 45
.2060	0.5	l	TNB	Twin-neck bottle	28 mm *
.2100	0.5	l	Sfs	Safety plastic bottle, black	DIN 45
.2110	0.5	l	Sfs	Safety plastic bottle black, with degassing lock	DIN 45
.2300	0.5	l	PFA	Perfluorokoxy bottle - reusable	33 mm *
.2301	0.5	l	PFA	Perfluorokoxy bottle - refilling	33 mm *
.2600	0.5	kg	Pfw	Plastic bottle, wide neck	DIN 60
.2700	0.5	l	Fl	Glass bottle	DIN 45
.2800	0.5	kg	Gl	Glass bottle, wide neck	DIN 54W
.3000	1.0	l	Sfl	Safety plastic bottle	DIN 45

Package key	Package size	Package unit	Package article	Package article in full	Thread
.3010	1.0	l	Sfl	Safety plastic bottle, with degassing lock	DIN 45
.3060	1,0	l	TNB	Twin-Neck-Bottle	28 mm *
.3100	1.0	l	Sfs	Safety plastic bottle, black	DIN 45
.3110	1.0	l	Sfs	Safety plastic bottle, black, with degassing lock	DIN 45
.3600	1.0	kg	Pfw	Plastic bottle, wide neck	DIN 80
.3700	1.0	l	Fl	Glass bottle	DIN 45
.3710	1.0	l	Fl	Glass bottle, with degassing lock	DIN 45
.3735	1.0	l	Fl	Glass bottle, lock with septum	DIN 45
.3755	1.0	l	Fl	Glass bottle, lock with septum, shrink-wrapped	DIN 45
.3800	1.0	kg	Gl	Glass bottle, wide neck	DIN 54W
.4000	2.5	l	Sfl	Safety plastic bottle	DIN 45
.4010	2.5	l	Sfl	Safety plastic bottle, with degassing lock	DIN 45
.4100	2.5	l	Sfs	Safety plastic bottle, black	DIN 45
.4110	2.5	l	Sfs	Safety plastic bottle, black, with degassing lock	DIN 45
.4600	2.5	kg	Pfw	Plastic bottle, wide neck	DIN 80
.4700	2.5	l	Fl	Glass bottle	DIN 45
.4710	2.5	l	Fl	Glass bottle, with degassing lock	DIN 45
.5000	5.0	l	Pka	Plastic canister	ov 51
.5010	5.0	l	Pka	Plastic canister, UN-approved, with degassing lock	ov 51
.5020	5.0	l	Pka	Plastic canister, UN-approved	ov 51
.5100	5.0	l	Sfs	Safety plastic bottle, black	DIN 45
.5110	5.0	l	Pkb	Plastic canister, blue, UN-approved, with degassing lock	ov 51
.5120	5.0	l	Pkb	Plastic canister, blue, UN-approved	ov 51
.5200	5.0	l	Pkk	Kautex canister, stackable	31 mm *
.5300	5.0	l	Sfl	Safety plastic bottle	DIN 45
.5310	5.0	l	Sfl	Safety plastic bottle, with degassing lock	DIN 45
.5350	5.0	l	FSfl	Safety plastic bottle, fluorinated	DIN 45
.5400	5.0	l	Lc	Labcube®	31,5 mm *
.5500	5.0	l	Blk	Metal canister	42 mm *
.5600	5.0	kg	Pfw	Bucket-yellow, with poly-flat bag insert	210 mm *

Package key	Package size	Package unit	Package article	Package article in full	Thread
.5900	5.0	l	Afl	Aluminium bottle	DIN 45
.6000	10.0	l	Pka	Plastic canister	ov 51
.6010	10.0	l	Pka	Plastic canister, UN-approved, with degassing lock	ov 51
.6020	10.0	l	Pka	Plastic canister, UN-approved	ov 51
.6110	10.0	l	Pkb	Plastic canister, blue, UN-approved, with degassing lock	ov 51
.6120	10.0	l	Pkb	Plastic canister, blue, UN-approved	ov 51
.6350	10.0	l	Blk	Metal canister	60 mm *
.6400	10.0	l	Lc	Labcube®	31 mm *
.6600	10.0	kg	Pfw	Bucket-yellow, with poly-flat bag insert	290 mm *
.6700	10.0	l	Pkl	Plastic canister, UN-approved, conductive	nv 61
.6800	various sizes		Tf	Dropper bottle	
.7000	20.0	l	Pka	Plastic canister	ov 61
.7010	20.0	l	Pka	Plastic canister, UN-approved, with degassing lock	ov 61
.7020	20.0	l	Pka	Plastic canister, UN-approved	ov 61
.7110	20.0	l	Pkb	Plastic canister, blue, UN-approved, with degassing lock	ov 61
.7120	20.0	l	Pkb	Plastic canister, blue, UN-approved	ov 61
.7400	20.0	l	Pc	Polycube	31 mm *
.7500	25.0	kg	Psk	Polyethylene bag in a cartonage	
.7600	25.0	kg	Ps	Polyethylene bag	
.7700	25.0	l	Glst	Glass bottle in polystyrene	DIN 70
.7900	25.0	l	Pka	Plastic canister, UN-approved	ov 61
.7910	25.0	l	Pka	Plastic canister, UN-approved, with degassing lock	ov 61
.7950	25.0	l	Pkb	Plastic canister, blue, UN-approved	ov 61
.7960	25.0	l	Pkb	Plastic canister, blue, UN-approved, with degassing lock	ov 61
.8010	30.0	l	Pka	Plastic canister, UN-approved, with degassing lock	ov 61
.8020	30.0	l	Pka	Plastic canister, UN-approved	ov 61
.8110	30.0	l	Pkb	Plastic canister, blue, UN-approved, with degassing lock	ov 61
.8120	30.0	l	Pkb	Plastic canister, blue, UN-approved	ov 61
.8800	600	l	IBC	Container	150 mm *
.8900	30.0	l	Ekm	Steel jug, reusable	2" & ¾" TriSure lock

Package key	Package size	Package unit	Package article	Package article in full	Thread
.8901	30.0	l	Ekm	Steel jug, refilling	2" & ¾" TriSure lock
.9300	60.0	l	Pka	Plastic canister, with lock QAA	ov 71
.9400	60.0	l	Pkb	Plastic canister blue, with lock QAA	ov 71
.9450	60.0	l	Pkl	Plastic canister, conductive	nv 71
.9600	500	l	Pt	Plastic tank on pipe pallet	150 mm *
.9700	1000	l	IBC	Container	150 mm *
.9800	100	l	Pf	Plastic barrel, UN-approved	S70 x 6 / S38 x 6
.9900	200	l	Pf	Plastic barrel, UN-approved	S70 x 6 / S56 x 4
.9950	200	l	Pf	Plastic barrel, UN-approved, conductive	S70 x 6 / S56 x 4
.9900	200	l	Pf	Fût en plastique avec certification UN	S70 x 6 / S56 x 4
.9950	200	l	Pf	Fût en plastique avec certification UN, conducteur	S70 x 6 / S56 x 4

* Size of filling opening